


# sempermed®

## supreme



**quick don™**

### Latex Sterile Surgical Gloves

- Surgical, 0.65 AQL
- Extra supple with great tactile sensitivity
- Fully anatomical with curved fingers and rolled cuff
- Excellent comfort and fit
- Micro-rough textured surface - superior instrument handling and grip
- Fast and easy to don: Quick Don™ inner coating
- Washed to minimize skin irritants
- Patented synthetic inner coating: US patent 6440499

REORDER INFO		Glove Color: Cream	
<b>SPFP550</b>	Size 5.5	50 pairs/box	6 boxes/case
<b>SPFP600</b>	Size 6	50 pairs/box	6 boxes/case
<b>SPFP650</b>	Size 6.5	50 pairs/box	6 boxes/case
<b>SPFP700</b>	Size 7	50 pairs/box	6 boxes/case
<b>SPFP750</b>	Size 7.5	50 pairs/box	6 boxes/case
<b>SPFP800</b>	Size 8	50 pairs/box	6 boxes/case
<b>SPFP850</b>	Size 8.5	50 pairs/box	6 boxes/case
<b>SPFP900</b>	Size 9	50 pairs/box	6 boxes/case

### RECOMMENDED USES



### PRODUCT SPECS



Rev08102021

