





SEMPERSHIELD® CR

Chloroprene EXAMINATION GLOVES



100 / 10 /

Examination Chloroprene Gloves

- **Thickness:** Competes with a 5-6 mil glove
- Fit similar to latex, strength similar to nitrile
- Excellent comfort - east to don
- Glove conforms to hand, reducing hand fatigue
- Excellent barrier and chemical resistance
- Superior grip in dry and wet environments
- Meets FDA 21 CFR Part 177 requirements for food handling

RECOMMENDED USES



PRODUCT SPECS



REORDER INFO		Glove Color: Jade	
SSCR102	X-Small	100 gloves/box	10 boxes/case
SSCR103	Small	100 gloves/box	10 boxes/case
SSCR104	Medium	100 gloves/box	10 boxes/case
SSCR105	Large	100 gloves/box	10 boxes/case
SSCR106	X-Large	100 gloves/box	10 boxes/case

SAMPLES: Available by request - ask your sales representative

For maximum protection, we recommend changing chloroprene gloves every 30 minutes



Rev02112022

